



LW-E8F

Scientific

Centrifuges

The 8-place LW Scientific E8 is an economical fixed speed centrifuge designed to fit any budget. This unit accommodates most test tubes from 3ml - 15ml. The wide stance eliminates wobble, while the suction-cupped feet prevent slipping. The brushless motor (110V only) is virtually maintenance-free and the lid port allows easy viewing and touching of specimens while in use.

Economical and compact!



Advantages

- Economical choice for students or small labs
- Accurate speeds for specific g-forces
- Proper separations without cell damage
- Quiet, vibration-free operation
- Space-saving design

Features

- Angled head rotor
- 8-tube capacity - 3ml to 15ml size tubes
- 3,300 rpm fixed speed - for blood only
- Auto-off 30-minute timer with bell
- Brushless motor (110V only) - maintenance free
- Blue LED power indicator light
- Suction-cupped feet prevent slipping
- 110V power (220V available)
- Maximum RCF=1,350 g

Includes

- Eight-place tube rotor
- Eight 15ml tube sleeves (TS-15ml)
- Eight 13x75mm tube sleeve inserts (TS-1375)

Warranty

One year on materials and workmanship

Dimensions

- Height: 10" (254mm)
- Length: 10" (254mm)
- Width: 10" (254mm)
- Weight: 11 lbs (5kg)

Model #	Description
LW-E8F	E8 Centrifuge, fixed speed