



Scientific

# M24

## Centrifuges

The LWS-M24 Centrifuge is a truly affordable microhematocrit centrifuge. This unit is designed for its ease of use in the laboratory. Spin hematocrits in less than 5 minutes at 12,000 rpm. The M24 has a microhematocrit reader disk provided with every unit. Read the tubes directly on the rotor without picking them up.

Spin hematocrits in less than 5 minutes at 12,000 rpm

### Advantages

- Fast results in just five minutes
- Spins 24 standard 75mm capillary tubes
- Ready-to-use...no add-on costs
- Quiet, vibration-free operation
- Reader disk included

### Design Features

- Impact-resistant safe-view lid
- All-metal body construction
- Cooling fan for motor compartment
- Aluminum rotor for low heat transfer
- Heavy duty high-speed motor
- Safety shut-off switch in lid
- Locking lid
- Wide mouth for easy clean up



### Specifications

- 12,000 rpm fixed speed
- Nominal RCF=14,800 g
- Brushed motor
- 110V AC (220V AC available)

- 24-Place Microhematocrit Centrifuge
- Affordable
- Results in Less Than 5 Minutes

### Controls

- Illuminated power switch
- Electric brake for quick stop
- Auto-off 15-minute timer with bell

### Includes

- 24-place hematocrit rotor
- Reader disk— read tubes directly on the rotor

### Warranty

One year on materials and workmanship

Model #	Description
LWS-M24	Microhematocrit Centrifuge

### Dimensions and Weight

- Height: 10" (235mm)
- Length: 13" (330mm)
- Width: 12" (304mm)
- Weight: 23 lbs. (10.4kg)