



# Pro Service Kit

## Pro Service Kit

Keep your optical equipment operating flawlessly with the new Pro Service Kit from LW Scientific. Who could know more about optical performance than a manufacturer of laboratory microscopes? Even a novice can perform expert scope service using the easy-to-understand picture manual in the Pro Service Kit. All tools are included for routine adjustments, alignments, and assembly, and the specialized cleaning supplies are gentle enough to protect the most delicate multi-coated optics. We know that a booklet of 25 lens papers wouldn't last 2 weeks in a lab, so we've included large-sized supplies that will last for months. Each component of the Pro Service Kit is available for individual sale, or choose the Standard Cleaning Kit for basic cleaning.



Standard Cleaning Kit Available as well

- A "Must Have" for Microscope Owners
- Specialized Tools for Adjusting and Aligning
- Detailed Picture Instruction Manual

### Pro Service Kit Includes:

- Air Duster - 8 oz.
- Lint-free Lens Paper - 280 sheets
- "No-residue" Lens Cleaning Fluid - 1 oz.
- Cotton-tipped Applicators - 100
- Microfiber Lens Cloth - 9x9 inch
- Soft-bristled Dusting Brush - 1
- Micro-Glide Gear Lubricant - 1 tube
- 3 Allen Wrenches - 1.5, 2.0, 2.5mm
- Friction Collar Wrench - double-sided
- Precision Screwdriver Set - 6 sizes
- Nylon Carry Case - 4x7x10 inches

### Standard Cleaning Kit Includes:

- Air Duster - 8 oz.
- Lint-free Lens Paper - 280 sheets
- "No-residue" Lens Cleaning Fluid - 1 oz.
- Cotton-tipped Applicators - 100

### Applications

- Compound Microscopes
- Stereo Microscopes
- Telescopes
- Cameras
- All glass and plastic lenses

Model #	Description
M-PSK	Pro Service Kit in Carry Bag
M-SCK	Standard Cleaning Kit

### Who Needs It?

- Every LW Scientific Distributor Sales Rep
- Biology Teachers and Homeschoolers
- Lab Technicians and Microscopists
- Astronomers
- Physicians and Veterinarians

### Instruction Manual

- Detailed Pictures for every tool use
- Compound and Stereoscope Service Techniques
- Technical Terms For Microscopy

#### Pro Service Kit - Dims:

Height: 10.0" (250mm) Width: 7.0" (177mm)  
Depth: 4" (103 mm) Weight: 1.9 lb. (0.9 kg)

#### Standard Cleaning Kit - Dims:

Height: 5.5" (140mm) Width: 8.75" (225mm)  
Depth: 5.5" (140 mm) Weight: 1.9 lb. (0.9 kg)