



Scientific

Turbo Mixers

Lab Equipment

LW Scientific now offers you two versions of the very popular Turbo Mixers. Ideal for the quick and complete mixing of samples, our cup and platform models will bring convenience and efficiency to your lab. They come with a variable speed control of up to 3000 RPM and can be activated to vibrate continuously or by pressure to the head. The high quality motor comes with a one year warranty and has been designed to offer quiet and vibration free vortex mixing action.

Characteristics

- Suction cup feet
- Metal construction
- Motor fan
- Compact and durable design
- Variable speed control
- Vibration-free
- High quality motor

Technical Specifications

Speed: 3000 RPM
110V / 60Hz (220V available)

Controls

On / off / touch switch
Power on light

Warranty

One year in parts and manufacture

Dimensions and Weight

Height: 170mm (6.7")
Width: 160mm (6.3")
Depth: 130mm (5.2")
Weight: 3.2Kg (7lbs)

LW-TM1000:
Cup head dimensions: 50 x 50 x 25mm
(2 x 2 x 1 inches)

LW-TM2000:
Platform head: 100mm diameter x 25mm height
(4" x 1")



Cup Turbo Mixer LW-TM1000

Also Available:
Platform Turbo Mixer LW-TM2000

- Vibration-free
- Variable Speed
- Continuous Operation or Touch Control

Model #	Description
LWTM1000	Turbo Mixer - Cup
LWTM2000	Turbo Mixer - Platform