



Scientific

Ultra

### Centrifuges

The 8-place LW Scientific Ultra has a unique, space-saving triangle shape to fit easily into a corner. This unit, made of sturdy aluminum alloy, accomodates most test tubes from 3ml-15ml. The wide stance eliminates wobble, while the suction-cupped feet prevent slipping. The brushless motor (110V only) is virtually maintenance-free and the clear lid allows easy viewing of specimens while in use. Choose from 4 different models to meet your laboratory separation needs.

4 Models to Meet Any Lab Need

#### Advantages

- Accurate speeds for specific g-forces
- Proper separations without cell damage
- Ideal for blood, urine and other fluids
- Quiet, vibration-free operation
- Unique space-saving triangle design

#### Features

- 8-tube capacity - 3ml to 15ml size tubes
- Auto-off 30-minute timer with bell
- Brushless motor (110V only) - maintenance free
- Blue LED power indicator light
- Suction-cupped feet prevent slipping
- All-metal one-piece body
- 110V power (220V available)

#### Ultra Tach Zero-Lock (LWU8TL)

- Digital tachometer - 4 digit display
- Variable Speed Control
- Lid locked in place until rotor comes to a complete stop.
- Speed display from 0-3,300 rpm
- Maximum RCF = 1350 g

#### Ultra Select (LWU8S)

- Selectable 4-position Speed Control
  - 3,300 rpm for blood - 1,350 g
  - 1800 rpm for urines and semen - 400 g
  - 1,500 rpm for fecals - 275 g
  - 800 rpm for other tests - 80 g
- Timer and lid safety shut-off switch

#### Ultra Variable (LWU8V)

- 3,300 rpm variable speed - blood, urines, etc.
- Timer and lid safety shut-off switch
- Maximum RCF=1,350 g



LWU8TL



LWU8V

- All-Metal One-Piece Design
- 4 Models to Meet Any Lab Need
- Locking Lid for Safety

#### Ultra Fixed (LWU8F)

- 3,300 rpm fixed speed - for blood only
- Timer and lid safety shut-off switch
- Maximum RCF=1,350 g

#### Includes

- Eight-place tube rotor
- Eight 15ml tube sleeves (TS-15ml)
- Eight 13x75mm tube sleeve inserts (TS-1375)

#### Warranty

One year on materials and workmanship

#### Dimensions

- Height: 9" (230mm)
- Length: 13" (330 mm)
- Width: 13" (330mm)
- Weight: 15.5 lbs (7kg)

Model #	Description
LWU8TL	Variable Tach 0-rpm Lock
LWU8S	4-speed Selectable
LWU8V	Variable Speed
LWU8F	Fixed Speed