

# Video Accessories

## 35mm Camera Adapters



Our 35mm camera adapters connect most of our microscopes with your camera and are compatible with all major brands of SLR cameras.

*Part# PHT25 Mount for use with Rev III  
PHT26 Mount for use with Observers  
PHT33 Mount for use with Ach/Par/  
Vision/M-Series/Z-2 Zoom*

## EZ Capture Card



To connect one of our video cameras to a desktop PC, just use the EZ Capture Card. It installs in minutes and allows you to capture both still and live video images.

*Part# BV-EZC*

## Digital Camera Adapters



Our digital camera adapter works with many of the most popular digital cameras on the market. Three different diameter adapter rings (28, 37, 52mm) are included.

*Part# PHT-DIG*

## TFT Screen



The 5.5" screen connects directly to all of our cameras providing a portable video screen. Set it on a table, or mount it to our boom stand.

*Part# TFT55*