



HemoPoint[®] H2

Accurate hemoglobin and hematocrit results with one simple test.

CLIA'88
WAIVED



STANBIO
LABORATORY

The World's Source For Clinical Diagnostics



Get precise hemoglobin
and hematocrit results
—both in the same test!

One test. Two results.

Fast.

Using a single drop of blood, Stanbio's **HemoPoint® H2** offers accurate results for both hemoglobin and hematocrit tests—in less than a minute! The easy-to-read touch screen displays the current time, battery status and allows for easy operation.

Advanced.

Utilizing proven microcuvette technology, the **HemoPoint® H2** features a proprietary “soft load” cuvette holder to minimize blood contamination of the meter. The hemoglobin result is flagged if it is outside the user-defined limits.

Easy to Use.

The **HemoPoint® H2**'s user-friendly sampling technique minimizes training time and its rechargeable battery makes it completely mobile. Recall the last 100 test results with the touch of a button and attach an optional printer to record patient's test results—immediately, or from stored memory!



HemoPoint® H2 Products

Catalog No. G3000-001

Photometer

Includes **HemoPoint® H2** Photometer, user guide, quick reference guide, optics cleaner, control cuvette, power adaptor, training CD and wall chart



Catalog No. 3010-100

100 Microcuvettes

Includes two containers of 50 cuvettes each

Catalog No. 3010-200

200 Microcuvettes

Includes four containers of 50 cuvettes each



Catalog No. 3060-601

Hgb Controls

Includes set of six dropper bottles:
2 x 1.5 mL-low, 2 x 1.5 mL-normal, 2 x 1.5 mL-high

Catalog No. 3050-001

Optics Cleaner

Catalog No. G3100-001

Printer, DPU-414

Thermal printer with power adaptor and cable

HemoPoint® H2 Specs

Methodology _____	Optical absorption photometry
Principle _____	Azidemethemoglobin
Calibration _____	Factory calibrated
Sample Type _____	Capillary, venous, or arterial blood
Sample Size _____	10 µL
Hematocrit Range _____	36 - 54% (calculated)
Linearity _____	0 - 23.5 g/dL
Within Run Precision CV _____	< 1%
Total Precision CV _____	< 1.2%



CPT CODE: 85018QW



STANBIO
LABORATORY

The World's Source For Clinical Diagnostics

- General & Special Chemistry
- Pregnancy
- Point-of-Care
- Rapid Tests
- Therapeutic Drug Monitoring
- Serology
- Occult Blood
- Urinalysis
- Clinical Instrumentation

Stanbio Laboratory
1261 North Main Street
Boerne, Texas 78006 USA

866.782.6246 (U.S. Toll Free)
Tel: 830.249.0772
Fax: 830.249.0851

email: stanbio@stanbio.com

www.stanbio.com

DISTRIBUTED BY